

Product Datasheet

NECAB3 Peptide - middle region (orb2009805)

Description	This is a synthetic peptide designed for use in combination with anti-NECAB3 Antibody(ARP57670_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	ESVEAQSRLLCGSRRAGRRALRSVSRSTWSPGSSDTGRSSEAEMQWRLQV
MW	44kDa
Uniprot ID	Q96P71
NCBI	NP_112509
Expiration Date	6 months from date of receipt.