

Product Datasheet

NDUFS8 Peptide - C-terminal region (orb2009811)

Description	This is a synthetic peptide designed for use in combination with anti-NDUFS8 Antibody (ARP61825_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	CKLCEAICPAQAITIEAEPRADGSRRTTRYDIDMTKCIYCGFCQEACPVD
MW	23kDa
Uniprot ID	O00217
NCBI	NP_002487
Expiration Date	6 months from date of receipt.