

Product Datasheet

NAT8L Peptide - C-terminal region (orb2009831)

Description	This is a synthetic peptide designed for use in combination with anti-NAT8L Antibody (ARP49151_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	VAAHKLYESLGFRHMGASDHVYVPGMTLSLAERLFFQVRYHRYRLQLREE
MW	15kDa
Uniprot ID	Q8N9F0
NCBI	NP_848652
Expiration Date	6 months from date of receipt.