

Product Datasheet

MAD2L1BP Peptide - N-terminal region (orb2010037)

Description	This is a synthetic peptide designed for use in combination with anti-MAD2L1BP Antibody (ARP61356_P0...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	SAAVPDLEWYEKSEETHASQIELLETSSTQEPLNASEAFCDRCMVPVVF
MW	30kDa
Uniprot ID	Q15013
NCBI	NP_055443
Expiration Date	6 months from date of receipt.