

Product Datasheet

Jph2 Peptide - N-terminal region (orb2010228)

Description	This is a synthetic peptide designed for use in combination with anti-Jph2 Antibody (ARP49509_P050),...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	QGKRHGLGIETKGRWLYKGEWTHGFKGRYGIRQSTNSGAKYEGTWNNGLQ
MW	75kDa
Uniprot ID	Q9ET78
NCBI	NP_067541
Expiration Date	6 months from date of receipt.