

Product Datasheet

IPPK Peptide - middle region (orb2010258)

Description	This is a synthetic peptide designed for use in combination with anti-IPPK Antibody(ARP57655_P050), ...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	KTLQVQMLDLLDIEGLYPLYNRVERYLEEFPEERKTLQIDGPYDEAFYQK
MW	56kDa
Uniprot ID	Q9H8X2
NCBI	NP_073592
Expiration Date	6 months from date of receipt.