

Product Datasheet

FBX032 Peptide - N-terminal region (orb2010738)

Description	This is a synthetic peptide designed for use in combination with anti-FBX032 Antibody (ARP52865_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	QQQLNNIQITRPAFKGLTFTDLPLCLQLNIMQR LSDGRDLVSLGQAAPDL
MW	25kDa
Uniprot ID	Q0VAQ6
NCBI	NP_680482
Expiration Date	6 months from date of receipt.