

## Product Datasheet

### EIF3F Peptide - N-terminal region (orb2010885)

|                            |                                                                                                                                                                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | EIF3F Peptide - N-terminal region                                                                                                                                                                                                                                       |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                                                                      |
| <b>Preservatives</b>       | Lyophilized powder                                                                                                                                                                                                                                                      |
| <b>Form/Appearance</b>     | Lyophilized powder                                                                                                                                                                                                                                                      |
| <b>Storage</b>             | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.                                                                                                                 |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                   |
| <b>Application notes</b>   | This is a synthetic peptide designed for use in combination with EIF3F Rabbit Polyclonal Antibody (orb585179). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings. |
| <b>Protein Sequence</b>    | VILASIVDSYERRNEGAARVIGTLLGTVDKHSVEVTNCFVPHNESEDEV                                                                                                                                                                                                                       |
| <b>MW</b>                  | 39kDa                                                                                                                                                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>O00303</b>                                                                                                                                                                                                                                                           |
| <b>NCBI</b>                | <b>NP_003745</b>                                                                                                                                                                                                                                                        |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                                                          |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)