

Product Datasheet

SH2D4A antibody (orb20109)

Description	Goat polyclonal antibody to SH2D4A
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SH2D4A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx. 55kDa band observed in Human Kidney lysates (calculated MW of 52.7kDa according to NP_071354.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HQKKESMKKKQDE
Clonality	Polyclonal
MW	52.7
Entrez	63898
Dilution Range	ELISA: 1:1000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb20109
(0.1ug/ml)
staining of
Human Ki...