

Product Datasheet

CTDSP2 Peptide - N-terminal region (orb2011095)

Description	This is a synthetic peptide designed for use in combination with anti-CTDSP2 Antibody(ARP52096_P050)...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	MEHGSITQARREDALVLTQGLVSKSSPKKPRGRNIFKALFCCFRAQHV
MW	31kDa
Uniprot ID	O14595
NCBI	NP_005721
Expiration Date	6 months from date of receipt.