

## Product Datasheet

### Coq5 Peptide - N-terminal region (orb2011153)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | This is a synthetic peptide designed for use in combination with anti-Coq5 Antibody (ARP48963_P050),...   |
| <b>Tested Applications</b> | WB  |
| <b>Preservatives</b>       | Lyophilized powder  |
| <b>Form/Appearance</b>     | Lyophilized powder  |
| <b>Storage</b>             | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Protein Sequence</b>    | RLSLQEKRAAETHFGFETVSEKEKGGKVYQVFQSVARKYDLMNDMMSLGI  |
| <b>MW</b>                  | 35kDa   |
| <b>Uniprot ID</b>          | <a href="#">Q4G064</a>  |
| <b>NCBI</b>                | <a href="#">NP_001034111</a>  |
| <b>Expiration Date</b>     | 6 months from date of receipt.  |