

Product Datasheet

APOL3 antibody (orb20114)

Description	Goat polyclonal antibody to APOL3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	APOL3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx. 40kDa band observed in lysates of cell lines DAUDI and MOLT4 (calculated MW of 44.3kDa according to NP_663615.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	ELTQIYQRLNPCHTH
Clonality	Polyclonal
MW	44.3; 36.5; 21.6
Entrez	80833
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20114
(1µg/ml)
staining of
DAUDI lysa...