

Product Datasheet

ANKRD54 Peptide - middle region (orb2011787)

Description	This is a synthetic peptide designed for use in combination with anti-ANKRD54 Antibody(ARP52644_P050...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	EVHALKRLRDSANANDVETVQQLLEDGADPCAADDKGRGTALHFASCNGND
MW	32kDa
Uniprot ID	Q6NXT1
NCBI	NP_620152
Expiration Date	6 months from date of receipt.