

## Product Datasheet

### FSCN3 Peptide - middle region (orb2012029)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-FSCN3 Antibody(ARP53741_P050),...
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	WGKFALNFCIELQGSNLLTVLAPNGFYMRADQSGTLLADSEITRECIWE
<b>MW</b>	55kDa
<b>Uniprot ID</b>	<b>Q9NQT6</b>
<b>NCBI</b>	<b>NP_065102</b>
<b>Expiration Date</b>	6 months from date of receipt.