

## Product Datasheet

### CFP Peptide - N-terminal region (orb2012704)

<b>Description</b>	This is a synthetic peptide designed for use in combination with anti-CFP antibody (Catalog #: ARP56...
<b>Tested Applications</b>	IHC, WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>MW</b>	48kDa
<b>Uniprot ID</b>	<a href="#">P27918</a>
<b>NCBI</b>	<a href="#">NP_002612</a>
<b>Expiration Date</b>	6 months from date of receipt.