

Product Datasheet

MPP3 Peptide - N-terminal region (orb2012724)

| | |
|----------------------------|---|
| Description | This is a synthetic peptide designed for use in combination with anti-MPP3 antibody (Catalog #: ARP5... |
| Tested Applications | WB |
| Preservatives | Lyophilized powder |
| Form/Appearance | Lyophilized powder |
| Storage | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| Note | For research use only |
| MW | 66kDa |
| Uniprot ID | Q13368 |
| NCBI | NP_001923 |
| Expiration Date | 6 months from date of receipt. |