

Product Datasheet

FAU Peptide (orb2013299)

Description	This is a synthetic peptide designed for use in combination with anti-FAU antibody (Catalog #: ARP54...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	14kDa
Uniprot ID	P62861
NCBI	NP_001988
Expiration Date	6 months from date of receipt.