

Product Datasheet

CBLN4 Peptide (orb2013850)

Description	This is a synthetic peptide designed for use in combination with anti-CBLN4 antibody (Catalog #: ARP...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	22kDa
Uniprot ID	Q9NTU7
NCBI	NP_542184
Expiration Date	6 months from date of receipt.