

Product Datasheet

SEPT-9 Peptide (orb2014027)

Description	This is a synthetic peptide designed for use in combination with anti-SEPT9(septin 9) antibody (Cata...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	37kDa
Uniprot ID	Q9UHD8
NCBI	NP_001106968
Expiration Date	6 months from date of receipt.