

Product Datasheet

GALNTL4 Peptide - middle region (orb2014690)

Description	This is a synthetic peptide designed for use in combination with anti-GALNTL4 antibody (Catalog #: A...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	69kDa
Uniprot ID	Q6P9A2
NCBI	NP_940918
Expiration Date	6 months from date of receipt.