


Product Datasheet

NDEL1 antibody (orb20147)

Description	Goat polyclonal antibody to NDEL1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NDEL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 38+39kDa bands observed in Human Brain (Cerebellum) lysates (calculated MW of 37.0kDa according to NP_001020750.1 and of 38.4kDa according to NP_110435.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	CRNFAKDQASRKS Y
Clonality	Polyclonal
MW	37.0; 38.4
Entrez	81565
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb20147
(0.1ug/ml)
staining of
Human Ce...