

Product Datasheet

ATP5IF1 Peptide - N-terminal region (orb2014916)

Description	This is a synthetic peptide designed for use in combination with anti-ATPIF1 antibody (Catalog #: AR...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	10kDa
Uniprot ID	Q9UII2
NCBI	NP_057395
Expiration Date	6 months from date of receipt.