

Product Datasheet

GDE1 Peptide (orb2015282)

Description	This is a synthetic peptide designed for use in combination with anti-GDE1 antibody (Catalog #: ARP4...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	38kDa
Uniprot ID	Q9NZC3
NCBI	NP_057725
Expiration Date	6 months from date of receipt.