

Product Datasheet

ZDHHC17 Peptide (orb2015329)

Description	This is a synthetic peptide designed for use in combination with anti-ZDHHC17 antibody (Catalog #: A...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	44kDa
Uniprot ID	Q8IUH5
NCBI	EAW97333
Expiration Date	6 months from date of receipt.