

Product Datasheet

PSMG1 antibody (orb20157)

Description	Goat polyclonal antibody to PSMG1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	DSCR2 / PSMG1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx. 30kDa band observed in lysates in Jurkat (calculated MW of 30.3kDa according to NP_982257.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KLMTTNEIQSNIYT
Clonality	Polyclonal
MW	32.9; 30.3; 32.1; 22.8
Entrez	8624
Dilution Range	ELISA: 1:64000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



orb20157
(0.3ug/ml)
staining of
Jurkat l...