

Product Datasheet

UBE2D2 Peptide - middle region (orb2016563)

Description	This is a synthetic peptide designed for use in combination with anti-EDG8 antibody (Catalog #: ARP4...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	16kDa
Uniprot ID	P62837
NCBI	NP_003330
Expiration Date	6 months from date of receipt.