

Product Datasheet

LIN7C Peptide - middle region (orb2016598)

Description	This is a synthetic peptide designed for use in combination with anti-LIN7C antibody (Catalog #: ARP...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	22kDa
Uniprot ID	Q9NUP9
NCBI	NP_060832
Expiration Date	6 months from date of receipt.