

Product Datasheet

ADARB1 Peptide - N-terminal region (orb2017715)

Description	ADARB1 Peptide - N-terminal region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	80kDa
Uniprot ID	P78563
NCBI	NP_001028221
Expiration Date	6 months from date of receipt.