

Product Datasheet

SIDT2 Peptide - N-terminal region (orb2017830)

Description	This is a synthetic peptide designed for use in combination with anti-SIDT2 antibody (Catalog #: ARP...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	94kDa
Uniprot ID	Q8NBJ9
NCBI	NP_001035545
Expiration Date	6 months from date of receipt.