

Product Datasheet

PTTG1 Peptide - N-terminal region (orb2018552)

Description	This is a synthetic peptide designed for use in combination with anti-PTTG1 antibody (Catalog #: ARP...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	22kDa
Uniprot ID	O95997
NCBI	NP_004210
Expiration Date	6 months from date of receipt.