

Product Datasheet

ZNF667 Peptide - middle region (orb2019306)

Description	This is a synthetic peptide designed for use in combination with anti-ZNF667 antibody (Catalog #: AR...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	70kDa
Uniprot ID	Q5HYK9
NCBI	NP_071386
Expiration Date	6 months from date of receipt.