

Product Datasheet

KCNA7 Peptide (orb2019585)

Description	This is a synthetic peptide designed for use in combination with anti-KCNA7 antibody (Catalog #: ARP...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	50kDa
Uniprot ID	Q96RP8
NCBI	NP_114092
Expiration Date	6 months from date of receipt.