

Product Datasheet

TAF2 Peptide - middle region (orb2020357)

Description	This is a synthetic peptide designed for use in combination with anti-TAF2 antibody (Catalog #: ARP3...
Tested Applications	ChIP, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	137kDa
Uniprot ID	Q6P1X5
NCBI	NP_003175
Expiration Date	6 months from date of receipt.