

Product Datasheet

ZC3H7B Peptide (orb2020750)

Description	This is a synthetic peptide designed for use in combination with anti-ZC3H7B antibody (Catalog #: AR...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	110kDa
Uniprot ID	Q9UGR2
NCBI	NP_060060
Expiration Date	6 months from date of receipt.