

Product Datasheet

MMP19 Peptide (orb2020852)

Description	This is a synthetic peptide designed for use in combination with anti-MMP19 antibody (Catalog #: ARP...
Tested Applications	IHC, WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	56kDa
Uniprot ID	Q99542
NCBI	NP_002420
Expiration Date	6 months from date of receipt.