

Product Datasheet

NATIVE Human C-REACTIVE Protein [Out of stock] (orb2022137)

Description	Purified Human C-reactive protein is an acute phase protein, produced by the liver that is used as a...
Tested Applications	ELISA
Preservatives	Liquid. Tris buffered saline, 5 mM CaCl ₂ .
Form/Appearance	Liquid. Tris buffered saline, 5 mM CaCl ₂ .
Storage	2°C to 8°C
Note	For research use only
Purity	> 99%
MW	114 kDa
Uniprot ID	P02741
NCBI	NP_000558.2
Expiration Date	6 months from date of receipt.