

## Product Datasheet

# NATIVE Human THYROGLOBULIN Protein [Out of stock] (orb2022258)

<b>Description</b>	NATIVE Human THYROGLOBULIN Protein [Out of stock]
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	Liquid. Phosphate buffered saline
<b>Form/Appearance</b>	Liquid. Phosphate buffered saline
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Purity</b>	SDS PAGE: > 96%
<b>MW</b>	660 kDa, confirmed by chromatography and electrophoresis. Under certain conditions the subunit molecular weight may be observed (330 kDa).
<b>Uniprot ID</b>	<a href="#">P01266</a>
<b>NCBI</b>	<a href="#">NP_003226.4</a>
<b>Expiration Date</b>	6 months from date of receipt.