

Product Datasheet

NATIVE Human PROSTATE SPECIFIC ANTIGEN Protein (orb2022265)

Description	NATIVE Human PROSTATE SPECIFIC ANTIGEN Protein
Tested Applications	ELISA
Preservatives	Liquid. Phosphate buffered saline
Form/Appearance	Liquid. Phosphate buffered saline
Storage	-20°C
Note	For research use only
Purity	> 40% by SDS-PAGE
MW	30 kD
Uniprot ID	P07288
NCBI	NP_001025218.1
Expiration Date	6 months from date of receipt.