

Product Datasheet

NATIVE Human PROSTATE SPECIFIC ANTIGEN Protein (orb2022266)

| | |
|----------------------------|---|
| Description | NATIVE Human PROSTATE SPECIFIC ANTIGEN Protein |
| Tested Applications | ELISA, WB |
| Preservatives | Liquid. Phosphate buffered saline |
| Form/Appearance | Liquid. Phosphate buffered saline |
| Storage | -20°C |
| Note | For research use only |
| Application notes | Application Info: Western Blotting: This protein may be used in Western blotting applications using suitable anti prostate specific antigen antibodies for detection and analysis. |
| Purity | > 96% by SDS-PAGE |
| MW | 30 kD |
| Uniprot ID | P07288 |
| NCBI | NP_001025218.1 |
| Expiration Date | 6 months from date of receipt. |