

Product Datasheet

NATIVE Human MYOGLOBIN Protein (orb2022273)

Description	NATIVE Human MYOGLOBIN Protein
Preservatives	Purified protein from human heart - liquid
Form/Appearance	Purified protein from human heart - liquid
Storage	Store at -20°C only.Storage in frost-free freezers is not recommended.This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Note	For research use only
Purity	SDS PAGE: > 99%
Uniprot ID	P02144
NCBI	NP_005359.1
Expiration Date	6 months from date of receipt.