

Product Datasheet

NATIVE MYCOPLASMA PNEUMONIAE Protein (orb2022274)

| | |
|----------------------------|--|
| Description | NATIVE MYCOPLASMA PNEUMONIAE Protein |
| Tested Applications | ELISA |
| Preservatives | Liquid. Phosphate buffered saline |
| Form/Appearance | Liquid. Phosphate buffered saline |
| Storage | -80°C |
| Note | For research use only |
| Application notes | Application Info: This product may be used in a variety of immunoassay formats or may be further purified to meet the requirements of a particular assay format. |
| Expiration Date | 6 months from date of receipt. |