

Product Datasheet

NATIVE Human LAMININ Protein [Out of stock] (orb2022277)

Description	NATIVE Human LAMININ Protein [Out of stock]
Tested Applications	ELISA, FA
Preservatives	TRIS buffered saline
Form/Appearance	TRIS buffered saline
Storage	-20°C
Note	For research use only
Purity	> 90% (SDS-PAGE)
MW	In electrophoresis a major band is observed at 60 kDa and a minor band at 170-190 kDa.
Expiration Date	6 months from date of receipt.