

Product Datasheet

NATIVE Human PLASMA KALLIKREIN Protein [Out of stock] (orb2022279)

Description	NATIVE Human PLASMA KALLIKREIN Protein [Out of stock]
Preservatives	Purified protein from human plasma - liquid
Form/Appearance	Purified protein from human plasma - liquid
Storage	Store at -20°C only.Storage in frost-free freezers is not recommended.This product should be stored undiluted.Avoid repeated freezing and thawing as this may denature the antibody.Should this product contain a precipitate we recommend microcentrifugation before use.
Note	For research use only
Purity	SDS PAGE: > 94%
Uniprot ID	P03952
NCBI	NP_000883.2
Expiration Date	6 months from date of receipt.