

## Product Datasheet

# NATIVE Human FIBRINOGEN Protein [Out of stock] (orb2022286)

<b>Description</b>	NATIVE Human FIBRINOGEN Protein [Out of stock]
<b>Tested Applications</b>	ELISA, FA
<b>Preservatives</b>	Liquid. Citrate buffered saline
<b>Form/Appearance</b>	Liquid. Citrate buffered saline
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Purity</b>	35-65% of protein, of which > 90% clottable protein. Virtually free from plasminogen and plasmin
<b>MW</b>	340 kDa
<b>Uniprot ID</b>	<a href="#">P02671</a>
<b>NCBI</b>	<a href="#">NP_000499.1</a>
<b>Expiration Date</b>	6 months from date of receipt.