

## Product Datasheet

# NATIVE Human FIBRIN DEGRADATION PRODUCT X Protein (orb2022287)

<b>Description</b>	NATIVE Human FIBRIN DEGRADATION PRODUCT X Protein
<b>Preservatives</b>	Purified protein from human serum - liquid
<b>Form/Appearance</b>	Purified protein from human serum - liquid
<b>Storage</b>	Store at -20°C only.Storage in frost-free freezers is not recommended.This product should be stored undiluted.Avoid repeated freezing and thawing as this may denature the antibody.Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Note</b>	For research use only
<b>Purity</b>	SDS PAGE: > 95%
<b>MW</b>	MW: 240 kD
<b>Uniprot ID</b>	<b>P02679</b>
<b>NCBI</b>	<b>NP_000500.2</b>
<b>Expiration Date</b>	6 months from date of receipt.