

Product Datasheet

NATIVE Human FIBRIN DEGRADATION PRODUCT E Protein (orb2022288)

Description	NATIVE Human FIBRIN DEGRADATION PRODUCT E Protein
Tested Applications	ELISA
Preservatives	Lyophilized. Glycine buffered saline
Form/Appearance	Lyophilized. Glycine buffered saline
Storage	-20°C
Note	For research use only
MW	MW: 50 kD
Uniprot ID	P02679
NCBI	NP_000500.2
Expiration Date	6 months from date of receipt.