

## Product Datasheet

# NATIVE Human FIBRIN DEGRADATION PRODUCT D-MONOMER Protein [Out of stock] (orb2022289)

<b>Description</b>	NATIVE Human FIBRIN DEGRADATION PRODUCT D- MONOMER Protein [Out of stock]
<b>Preservatives</b>	Purified protein from human serum - liquid
<b>Form/Appearance</b>	Purified protein from human serum - liquid
<b>Storage</b>	Store at +4°C or at -20°C if preferred.Storage in frost-free freezers is not recommended.This product should be stored undiluted.
<b>Note</b>	For research use only
<b>Purity</b>	SDS PAGE: > 95%
<b>MW</b>	45 kDa under reducing conditions and 94 kDa under non reducing conditions.
<b>Uniprot ID</b>	<a href="#">P02671</a>
<b>NCBI</b>	<a href="#">NP_000499.1</a>
<b>Expiration Date</b>	6 months from date of receipt.