

Product Datasheet

NATIVE Human FACTOR H Protein (orb2022292)

Description	NATIVE Human FACTOR H Protein
Tested Applications	IHC-Fr, WB
Immunogen	Purified guinea pig calretinin
Preservatives	Liquid. TRIS buffered saline
Form/Appearance	Liquid. TRIS buffered saline
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	SDS PAGE: > 95%
MW	150 kDa
Uniprot ID	P08603
NCBI	NP_000177.2
Expiration Date	6 months from date of receipt.