

Product Datasheet

NATIVE CARBONIC ANHYDRASE Protein (orb2022294)

Description	NATIVE CARBONIC ANHYDRASE Protein
Tested Applications	ELISA
Immunogen	E. coli strain CVL Moon 987P O:9K103.
Preservatives	Liquid. Phosphate buffered saline
Form/Appearance	Liquid. Phosphate buffered saline
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P00921
NCBI	NP_848667.1
Expiration Date	6 months from date of receipt.