

## Product Datasheet

# NATIVE Human C-REACTIVE Protein [Out of stock] (orb2022297)

<b>Description</b>	NATIVE Human C-REACTIVE Protein [Out of stock]
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	Liquid. TRIS buffered saline
<b>Form/Appearance</b>	Liquid. TRIS buffered saline
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Purity</b>	SDS PAGE: > 96%
<b>Source</b>	Human Plasma
<b>MW</b>	114 kDa
<b>Uniprot ID</b>	<b>P02741</b>
<b>NCBI</b>	<b>NP_000558.2</b>
<b>Expiration Date</b>	6 months from date of receipt.